

Chemical Exchange Saturation Transfer MRI

Contrast and Quantification

Paula Croal, PhD
paula.croal@eng.ox.ac.uk

*Institute of Biomedical Engineering & Wellcome Centre for Integrative Neuroimaging,
University of Oxford, Oxford, United Kingdom*




CEST MRI: The Basics

- Chemical Exchange Saturation Transfer MRI
- Molecular imaging technique
- Detect low concentration of metabolites that have exchangeable protons
 - Amides
 - Amines
 - Hydroxyls
 - Creatine


Why CEST?

- Means of enhancing the sensitivity of MRI to broad range of solute molecules.
- Non-invasive means to measure biophysical and physiological properties.
- Clinically feasible and increasing evidence of clinical role.
- Broad range of CEST MRI depending on application.
 - Endogenous (e.g. APT CEST)
 - Exogenous (e.g. GlucoseCEST)


A Typical MRI Experiment




A Typical MRI Experiment: Saturation



The CEST Effect





The CEST Effect: Exchange



APT CEST MRI


What are we trying to measure?

- Amide Proton Transfer (APT) CEST MRI
- In simplest terms, we just want to know what the effect size is at 3.5 ppm in relation to water
- A key appeal of the amide pool is that the exchange rate (k_{amide}) is base-catalysed
 - pH measurement
- Broad range of clinical applications






Measuring the APT Effect

- Asymmetry Analysis



$$MTR_{asym}(\Delta) = \frac{S(-\Delta) - S(\Delta)}{S_0}$$



Ma et al., JMRI, 2016


Zhao et al, MRM, 2012

Wang et al., Chin Med J, 2015

Croal et al., ISMRM, 2019



Is it really that simple?

- The Z-spectra is actually made up from a range of signal contributions
- Direct water saturation
- Multiple pools
- Macromolecular Effects (non-symmetrical)
- B_0 offsets
- T_1/T_2 effects
- Acquisition parameters





Correcting for B_0 Effects

- Magnetic field inhomogeneities can result in a shifting of the spectra
- Rather problematic if left uncorrected
- Solution is to estimate your field shift and sample spectrum accordingly




Macromolecular effects

- You can ignore the fact that they are asymmetrical...
 - Proceed with asymmetry analysis
- You could try and incorporate them into your analysis.
 - Both model-free and model-based options




Accounting for Multiple Pools

- The advantage of incorporating macromolecular effects into your analysis, is that you can also use this approach to account for additional pools.
- The challenge is how best to characterise pool(s) of interest, while accounting for confounding pools.
- Your pool of interest will depend on what physiological parameter you are trying to measure.



Lorentzian Difference Analysis



Overview:

- Z-spectrum described assumed to be:
 - Water pool (Lorentzian)
 - Solute pools of interest
- Exploratory
- Doesn't require a whole z-spectra.


Challenges:

- MT affects need to be well-accounted for.
- No T_1 correction.
- Non-quantitative in physiological terms.

Lorentzian Difference Analysis: Tools


Quantiphyse

- Visualisation and quantification software
- Freely available for research purposes
- www.quantiphyse.org




Bayesian Model-based Analysis

- The Bloch-McConnell equations allow us to characterise both the MR signal characteristics of water, as well as the exchange of protons between pools of interest and water.
- The user has control over the number of pools modelled; initialising parameters such as T_1 , T_2 , proton concentration (M_0), and exchange rate (k) to literature or experimental values.




Utilising a Bayesian model allows us to incorporate image priors, and attach degree of certainty

Bayesian Model-based Analysis




$$CESTR^* = \frac{S_{Water}(\Delta\omega) - S_{water+pool\ of\ interest}(\Delta\omega)}{M_0}$$



Inherent
B₀
correction


Bayesian Model-based Analysis




Bayesian Model-Based Analysis: CESTR*

$$CESTR^* = \frac{S_{Water}(\Delta\omega) - S_{water+pool\ of\ interest}(\Delta\omega)}{M_0}$$

- Model-based measures of exchange rate and concentration are not independent.
- Current approach is to combine into a ratio metric
 - Serves as an internal reference point
 - Comparable across subjects
 - Quantitative, yet comparable to MTR_{asym}
- Shown to be more robust than “conventional” metrics
- However, a solution would be to acquire at multiple powers



Bayesian Model-based Analysis



Overview:

- Pools initialised to literature and/or experimental values
- BM Equations used to model full Z-spectra
- ‘Pure’ CEST effect than characterised by CESTR*


Challenges:

- Relatively large time requirement
 - Acquisition
 - Analysis
- Challenging to implement (**Quantiphyse**)
- Reliance on metrics?

Bayesian Model-Based Analysis : Tools



Quantiphyse

- Visualisation and quantification software
- Freely available for research purposes
- www.quantiphyse.org




Does Analysis Method Really Matter?

- There are a vast array of CEST metrics available to use.
- Your choice will likely be influenced by a variety of factors.
 - Experimental design
 - Clinical population
 - Availability of resources
- However, choice of metric may significantly impact interpretation.



Is it Really That Simple?

Revisiting Asymmetry Analysis



$$MTR_{asym}(\Delta) = \frac{S(-\Delta) - S(\Delta)}{S_0}$$

Overview:

- Aims to isolate CEST contrast.
- Only requires 3 frequency offsets.
- Straightforward analysis.
- Broad range of applications.

Challenges:

- Pools on both sides.
- Asymmetric MT effects.
- No T_1 correction.

Conclusions:

- As with many imaging techniques, there is no consensus on how best to quantify parameters of interest.
- General trend to use “metrics”, however there are a large number available.
 - Choice of metric may impact interpretation
- Model-based analysis likely offers a robust approach.
 - **Bayesian model-based (CESTR*)**
 - Inherent B0 correction
 - Ability to incorporate “image priors”
 - May be particularly important in pathology

Online Materials

<https://quantiphyse.readthedocs.io/en/latest/cest/cest.html>

- General user guide
- Tutorials: Preclinical and clinical

Practical Session: Quantiphyse

- Model-based analysis of APT CEST MRI
- Choice of preclinical ischemia or clinical brain tumour
- Run from virtual machine

Acknowledgements:

Quantiphyse Team

Michael Chappell

Martin Craig

Benjamin Irving

Moss Zhao

Alex Smith

Flora Kennedy-McConnell

Yunus Mysaib

Tom Kirk

Nicholas Blockley

Matthew Cherukara

Prof. Sir. Mike Brady

Prof. Julia Schnabel



CANCER
RESEARCH
UK

EPSRC
Engineering and Physical Sciences
Research Council

CANCER IMAGING CENTRE

SUPPORTED BY

NIHR | National Institute
for Health Research